

## TUMOR CELLS



### ISOLATION OF TUMOR CELLS WITH COLLAGENASE NB

In the following, application details for isolation of cells from tumor tissue from breast, uterus, thyroid gland, lung, and colon are presented.

#### Tumor tissue from breast, uterus, thyroid gland or lung:

Species	Enzyme	Conditions	
Human*	Collagenase NB 8	3.5-5.5 PZ U/ml, 3-5 h, 37 °C	

#### Tumor tissue from colon:

Species	Enzyme	Conditions	
Human*	Collagenase NB 4G	0.5-1 PZ U/ml, RT	

\* For information on GMP grade Collagenase NB please [contact us!](#)

**Collagenase NB 8 Broad Range (cat. no. S1745601)** is characterized by a high collagenase activity ( $\geq 0.9$  PZ U/mg) and reduced levels of proteolytic side activities.

**Collagenase NB 4G Proved Grade (cat. no. S1746501)** has a collagenase activity of  $\geq 0.18$  PZ U/mg and a balanced mix of proteolytic side activities.

#### Isolation conditions

The suggested Collagenase NB concentrations and isolation conditions should be regarded as starting points. Collagenase NB dosage is recommended in PZ units to ensure consistent results.

For further information on optimization of tissue digestion please contact us at (631) 348-0333 or [collagenase@crescentchemical.com](mailto:collagenase@crescentchemical.com) or visit us at [www.crescentchemical.com](http://www.crescentchemical.com).

#### Ordering:

Crescent Chemical Co., Inc.  
2 Oval Dr.  
Islandia, N.Y. 11749  
Phone: (631)348-0333 Fax: (631) 348-0913  
E-Mail: [collagenase@crescentchemical.com](mailto:collagenase@crescentchemical.com)  
[www.crescentchemical.com](http://www.crescentchemical.com)